

Anti-LPCAT1 [R07-4K3] Monoclonal Antibody

AGMB00523

Description

This Anti-LPCAT1 [R07-4K3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00523
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	AYTL2, lpcat, AGPAT9, PFAAP3, AGPAT10, LPCAT-1, lysoPAFAT
Applications:	WB IHC-F IHC-P ICC/IF FC IP
Research Area:	-
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-4K3
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8NF37
Immunogen:	Recombinant protein of human LPCAT1
Gene ID:	79888
Gene Name:	LPCAT1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 59 kDa, Observed MW: 59 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:200-1:2000
	ICC/IF	1:100-1:200
	FC	1:20-1:50
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.