

Anti-SIRT6 [R04-9G5] Monoclonal Antibody

AGMB00536

Description

This Anti-SIRT6 [R04-9G5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00536
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	2810449N18Rik, AI043036, Mono ADP ribosyltransferase sirtuin 6, NAD-dependent protein deacetylase sirtuin-6, Regulatory protein SIR2 homolog 6, Regulatory protein SIR2 homolog, SIR2 like 6, SIR2 like protein 6, Sir2 related protein type 6, SIR2-like protein 6, SIR2
Applications:	WB ICC/IF IP
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R04-9G5
Clonality:	Monoclonal Antibody
SwissProt ID:	Q8N6T7
Immunogen:	A synthetic peptide of human SIRT6

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	51548
Gene Name:	SIRT6
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 39 kDa, Observed MW: 39 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	ICC/IF	1:500-1:1000
	IP	1:20-1:50
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	