

Anti-LYRIC/AEG1 [R02-7D4] Monoclonal Antibody

AGMB00575

Description

This Anti-LYRIC/AEG1 [R02-7D4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00575
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MTDH, AEG1, LYRIC, Protein LYRIC, 3D3/LYRIC, Astrocyte elevated gene-1 protein, AEG-1, Lysine-rich CEACAM1 co-isolated protein, Metadherin, Metastasis adhesion protein, 3D3, 3D3/LYRIC, AEG 1,
Applications:	WB IHC-F IHC-P ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R02-7D4
Clonality:	Monoclonal Antibody
SwissProt ID:	Q86UE4
Immunogen:	A synthetic peptide of human LYRIC/AEG1
Gene ID:	92140

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: MTDH
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 64 kDa, Observed MW: 75 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:2000-1:5000
	IHC-P	1:100-1:200
	ICC/IF	1:200-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.