

Anti-F4/80 [R08-4K-9] Monoclonal Antibody

AGMB00608

Description

This Anti-F4/80 [R08-4K-9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00608
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Adhesion G protein-coupled receptor E1, Adgre1, Cell surface glycoprotein F4/80, EGF-like module receptor 1, Adgre1, Emr1, Gpf480
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R08-4K-9
Clonality:	Monoclonal Antibody
SwissProt ID:	Q61549
Immunogen:	Peptide
Gene ID:	13733

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: Adgre1

Host Species: Rabbit

Isotype: IgG

Purification: Protein A

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.