

Anti-WIP1 [R04-118] Monoclonal Antibody

AGMB00622

Description

This Anti-WIP1 [R04-118] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00622
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	WIP1-1, WIP1 49 kDa, WIP1 1, WD repeat domain phosphoinositide-interacting protein 1, Atg18 protein homolog
Applications:	WB IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R04-118
Clonality:	Monoclonal Antibody
SwissProt ID:	Q5MNZ9
Immunogen:	A synthesized peptide derived from human WIP1
Gene ID:	55062

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: WIPI1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 49 kDa, Observed MW: 49 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	