

Anti-Bcl2 A1 [R09-7G1] Monoclonal Antibody

AGMB00672

Description

This Anti-Bcl2 A1 [R09-7G1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00672
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	BCL2A1, ACC-1, ACC-2, Bcl-2-like protein 5, Bcl-2-related protein A1, BCL2-related protein A1, BFL1, Bcl2-L-5, BCL2L5, GRS, HBPA1, Protein BFL-1, Protein GRS
Applications:	WB IHC-P FC
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R09-7G1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q16548
Immunogen:	A synthesized peptide derived from human BCL2A1
Gene ID:	597

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: BCL2A1

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 20 kDa, Observed MW: 27 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.