

Anti-TRIM25 [R06-1H2] Monoclonal Antibody

AGMB00783

Description

This Anti-TRIM25 [R06-1H2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00783
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	TRIM25, EFP, RNF147, ZNF147, E3 ubiquitin/ISG15 ligase TRIM25, Estrogen-responsive finger protein, RING finger protein 147, Tripartite motif-containing protein 25, Ubiquitin/ISG15-conjugating enzyme TRIM25, Zinc finger protein 147
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-1H2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q14258
Immunogen:	A synthetic peptide of human TRIM25

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	7706
Gene Name:	TRIM25
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 71 kDa, Observed MW: 71 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:1000-1:5000
IHC-P	1:100-1:200
ICC/IF	1:200-1:1000
FC	1:50-1:100
IP	1:20-1:100

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.