

# Anti-E2F2 [R03-6G9] Monoclonal Antibody

## AGMB00788

### Description

---

This Anti-E2F2 [R03-6G9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                                 |
|-----------------------|-------------------------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB00788                                                                                       |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                             |
| <b>Synonyms:</b>      | De2f2, E2F transcription factor 2, E2F-2, Transcription factor E2F2, E2F transcription factor 2 |
| <b>Applications:</b>  |              |
| <b>Research Area:</b> | Epigenetics                                                                                     |
| <b>Form:</b>          | Liquid                                                                                          |

### Antibody Data

---

|                      |                                               |
|----------------------|-----------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                             |
| <b>Clone:</b>        | R03-6G9                                       |
| <b>Clonality:</b>    | Monoclonal Antibody                           |
| <b>SwissProt ID:</b> | Q14209                                        |
| <b>Immunogen:</b>    | A synthesized peptide derived from human E2F2 |
| <b>Gene ID:</b>      | 1870                                          |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Gene Name:</b>        | E2F2                                       |
| <b>Host Species:</b>     | Rabbit                                     |
| <b>Isotype:</b>          | IgG                                        |
| <b>Purification:</b>     | Affinity Chromatography                    |
| <b>Conjugated:</b>       | Unconjugated                               |
| <b>Modification:</b>     | Unmodified                                 |
| <b>Molecular Weight:</b> | Calculated MW: 48 kDa, Observed MW: 48 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| <b>Antibody Dilution Ratio:</b>           | Application                                                                                                                                                                                                                                             | Antibody Dilution Ratio |
|-------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------|
|                                           | WB                                                                                                                                                                                                                                                      | 1:1000-1:2000           |
|                                           | FC                                                                                                                                                                                                                                                      | 1:50-1:100              |
| <b>Protein Quantification (Optional):</b> | To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <a href="https://www.assaygenie.com/bradford-protein-assay-protocol/">https://www.assaygenie.com/bradford-protein-assay-protocol/</a> to view the full protocol. |                         |