

# Anti-Desmoglein 2/DSG2 [R03-5P2] Monoclonal Antibody

## AGMB00798

### Description

---

This Anti-Desmoglein 2/DSG2 [R03-5P2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                            |
|-----------------------|----------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB00798                                                                  |
| <b>Contents:</b>      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                        |
| <b>Synonyms:</b>      | desmoglein-2, DSG2, HDGC, ARVD10, CDHF5, CMD1BB, ARVC 10, ARVC10, ARVD 10, |
| <b>Applications:</b>  | <b>WB</b> <b>IHC-P</b> <b>ICC/IF</b> <b>FC</b>                             |
| <b>Research Area:</b> | Signal Transduction                                                        |
| <b>Form:</b>          | Liquid                                                                     |

### Antibody Data

---

|                      |                                                       |
|----------------------|-------------------------------------------------------|
| <b>Reactivity:</b>   | Human, Mouse, Rat                                     |
| <b>Clone:</b>        | R03-5P2                                               |
| <b>Clonality:</b>    | Monoclonal Antibody                                   |
| <b>SwissProt ID:</b> | Q14126                                                |
| <b>Immunogen:</b>    | A synthesized peptide derived from human Desmoglein 2 |
| <b>Gene ID:</b>      | 1829                                                  |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** DSG2

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 122 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | WB                 | 1:1000-1:2000                  |
|                                 | IHC                | 1:100-1:200                    |
|                                 | IF                 | 1:50-1:200                     |
|                                 | FC                 | 1:50-1:100                     |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.