

Anti-RUNX2 [R08-7F3] Monoclonal Antibody

AGMB00805

Description

This Anti-RUNX2 [R08-7F3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00805
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CCD, AML3, CCD1, CLCD, OSF2, CBFA1, OSF-2, PEA2aA, PEBP2aA, CBF-alpha-1
Applications:	IHC-P ICC/IF
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R08-7F3
Clonality:	Monoclonal Antibody
SwissProt ID:	Q13950
Immunogen:	A synthesized peptide derived from human RUNX2
Gene ID:	860

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: RUNX2

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 56 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	