

Anti-FADD [R02-1A6] Monoclonal Antibody

AGMB00880

Description

This Anti-FADD [R02-1A6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00880
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FADD, MORT1, GIG3, Protein FADD, FAS-associated death domain protein, FAS-associating death domain-containing protein, Growth-inhibiting gene 3 protein, Mediator of receptor induced toxicity, FADD protein, FADD_HUMAN,
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R02-1A6
Clonality:	Monoclonal Antibody
SwissProt ID:	Q13158
Immunogen:	A synthetic peptide of human FADD

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	8772
Gene Name:	FADD
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 23 kDa, Observed MW: 28 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:1000-1:5000
IHC	1:20-1:100
ICC/IF	1:20-1:50
FC	1:20
IP	1:20-1:50

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.