

Anti-MCL1 [R03-3Y-5] Monoclonal Antibody

AGMB00937

Description

This Anti-MCL1 [R03-3Y-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00937
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MCL1, BCL2L3, Induced myeloid leukemia cell differentiation protein Mcl-1, Bcl-2-like protein 3, Bcl2-L-3, Bcl-2-related protein EAT/mcl1, mcl1/EAT, Bcl 2 related protein EAT/mcl1, Bcl-2-like protein 3, Bcl-2-related protein EAT/mcl1,
Applications:	IHC-P
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-3Y-5
Clonality:	Monoclonal Antibody
SwissProt ID:	Q07820
Immunogen:	A synthesized peptide derived from human MCL1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 4170
Gene Name: MCL1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.