

Anti-CD69 [R02-5E-4] Monoclonal Antibody

AGMB00944

Description

This Anti-CD69 [R02-5E-4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB00944
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD69, CLEC2C, Early activation antigen CD69, Activation inducer molecule, AIM, BL-AC/P26, C-type lectin domain family 2 member C, EA1, Early T-cell activation antigen p60, GP32/28, Leukocyte surface antigen Leu-23, MLR-3, CD antigen CD69, Activation inducer molecule (AIM/CD69), Activation inducer molecule,
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R02-5E-4
Clonality:	Monoclonal Antibody
SwissProt ID:	Q07108
Immunogen:	Recombinant protein of human CD69

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 969
Gene Name: CD69
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.