

Anti-CYP27A1 [R05-9I3] Monoclonal Antibody

AGMB01002

Description

This Anti-CYP27A1 [R05-9I3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01002
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CYP27A1, CYP27, Sterol 26-hydroxylase, mitochondrial, 5-beta-cholestane-3-alpha, 7-alpha, 12-alpha-triol 27-hydroxylase, Cytochrome P-450C27/25, Cytochrome P450 27, Sterol 27-hydroxylase, Vitamin D(3) 25-hydroxylase, 12-alpha-triol 27-hydroxylase, 5-beta-cholestane-3-alpha, 5-beta-cholestane-3-alpha, 7-alpha, 12-alpha-triol 26-hydroxylase,
Applications:	WB IHC-P ICC/IF FC
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R05-9I3
Clonality:	Monoclonal Antibody
SwissProt ID:	Q02318
Immunogen:	A synthesized peptide derived from human CYP27A1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 1593

Gene Name: CYP27A1

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 60 kDa, Observed MW: 50 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.