

# Anti-hnRNP U/p120 [R05-7N7] Monoclonal Antibody

## AGMB01028

### Description

---

This Anti-hnRNP U/p120 [R05-7N7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB01028
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	Heterogeneous nuclear ribonucleoprotein U, hnRNP U protein, HNRNPU, hnRNPU protein, HNRPU, HNRPU_HUMAN, p120, p120 nuclear protein, SAF A, SAF-A, SAFA, Scaffold attachment factor A, U21.1., Heterogeneous nuclear ribonucleoprotein U, hnRNP U, hnRNP U protein,
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b> <b>ICC/IF</b> <b>FC</b>
<b>Research Area:</b>	Epigenetics
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	R05-7N7
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q00839
<b>Immunogen:</b>	A synthesized peptide derived from human hnRNP U

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 3192

**Gene Name:** HNRNPU

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 91 kDa, Observed MW: 120 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

**Antibody Dilution Ratio:**

Application	Antibody Dilution Ratio
WB	1:1000-1:2000
IHC	1:100-1:200
IF	1:50-1:200
FC	1:50-1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.