

Anti-TFIIB [R06-4C2] Monoclonal Antibody

AGMB01036

Description

This Anti-TFIIB [R06-4C2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01036
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Transcription Initiation Factor IIB, General transcription factor IIB, General transcription factor TFIIB, gtf2b, RNA polymerase II transcription factor IIB, S300 II, S300-II, TF IIB, TF2B, TF2B_HUMAN, TFIIB,
Applications:	WB
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-4C2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q00403
Immunogen:	A synthetic peptide of human TFIIB
Gene ID:	2959

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: GTF2B
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 35 kDa, Observed MW: 35 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:1000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.