

Anti-SUMO Conjugating Enzyme UBC9 [R01-6I6] Monoclonal Antibody

AGMB01093

Description

This Anti-SUMO Conjugating Enzyme UBC9 [R01-6I6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01093
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	UBE2I, UBC9, UBCE9, SUMO-conjugating enzyme UBC9, SUMO-protein ligase, Ubiquitin carrier protein 9, Ubiquitin carrier protein I, Ubiquitin-conjugating enzyme E2 I, Ubiquitin-protein ligase I, p18
Applications:	WB IHC-F IHC-P ICC/IF FC IP
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-6I6
Clonality:	Monoclonal Antibody
SwissProt ID:	P63279
Immunogen:	A synthetic peptide of human UBE2I/UBC9
Gene ID:	7329

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	UBE2I
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 18 kDa, Observed MW: 18 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:1000-1:5000
IHC-P	1:100-1:200
ICC/IF	1:50-1:100
FC	1:20-1:50
IP	1:20

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.