

Anti-RHOA/B/C [R66-2G-8] Monoclonal Antibody

AGMB01145

Description

This Anti-RHOA/B/C [R66-2G-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01145
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ARH12, ARHA, H12, RHO12, RHOA, RHOB, RHOC, RHOH12, RHOC, Rho-related GTP-binding protein RhoC, Small GTP binding protein RhoA, Transforming protein RhoA
Applications:	WB IHC-P ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R66-2G-8
Clonality:	Monoclonal Antibody
SwissProt ID:	P61586/P62745/P08134
Immunogen:	A synthesized peptide derived from human Rho
Gene ID:	387

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: RHOA

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 22 kDa, Observed MW: 22 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.