

Anti-RAP1A+RAP1B [R03-1K4] Monoclonal Antibody

AGMB01151

Description

This Anti-RAP1A+RAP1B [R03-1K4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01151
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	C21KG, G 22K, K REV, KREV 1, KREV, KREV1, RAL1B, RAP 1, RAP 1A, RAP 1B, RAP1, RAP1B, SMGP21
Applications:	WB FC IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R03-1K4
Clonality:	Monoclonal Antibody
SwissProt ID:	P61224/P62834
Immunogen:	A synthesized peptide derived from human RAP1A + RAP1B
Gene ID:	5908/5906

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: RAP1B/RAP1A
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 21 kDa, Observed MW: 21 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	FC	1:50-1:100
	IP	1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.