

Anti-BAP31 [R03-9F1] Monoclonal Antibody

AGMB01277

Description

This Anti-BAP31 [R03-9F1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01277
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	BCAP31, BAP31, DXS1357E, B-cell receptor-associated protein 31, BCR-associated protein 31, Bap31, 6C6-AG tumor-associated antigen, Protein CDM, p28
Applications:	WB IHC-F IHC-P ICC/IF
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	R03-9F1
Clonality:	Monoclonal Antibody
SwissProt ID:	P51572
Immunogen:	A synthetic peptide of human BAP31
Gene ID:	10134

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: BCAP31

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 28 kDa, Observed MW: 28 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:200-1:500
	ICC/IF	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.