

Anti-SMARCA2 / BRM [R02-3E2] Monoclonal Antibody

AGMB01280

Description

This Anti-SMARCA2 / BRM [R02-3E2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01280
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	SMARCA2, BAF190B, BRM, SNF2A, SNF2L2, Probable global transcription activator SNF2L2, ATP-dependent helicase SMARCA2, BRG1-associated factor 190B, BAF190B, Protein brahma homolog, hBRM, SNF2-alpha, SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily A member 2
Applications:	WB IHC-P ICC/IF FC
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse
Clone:	R02-3E2
Clonality:	Monoclonal Antibody
SwissProt ID:	P51531
Immunogen:	A synthetic peptide of human SMARCA2/BRM

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 6595

Gene Name: SMARCA2

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 181 kDa, Observed MW: 190 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:200-1:2000
	ICC/IF	1:20-1:50
	FC	1:200-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.