

Anti-Nampt [R05-4C4] Monoclonal Antibody

AGMB01399

Description

This Anti-Nampt [R05-4C4] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01399
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Nicotinamide phosphoribosyltransferase, NAMPRtase, Nampt, Pre-B-cell colony-enhancing factor 1, Pre-B cell-enhancing factor, Visfatin
Applications:	WB IHC-F IHC-P ICC/IF FC IP
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R05-4C4
Clonality:	Monoclonal Antibody
SwissProt ID:	P43490
Immunogen:	A synthetic peptide of human Visfatin
Gene ID:	10135

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NAMPT
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 56 kDa, Observed MW: 56 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:200-1:2000
	ICC/IF	1:100
	FC	1:200-1:500
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.