

Anti-CD80 [R03-2J8] Monoclonal Antibody

AGMB01541

Description

This Anti-CD80 [R03-2J8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01541
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD80, CD28LG, CD28LG1, LAB7, T-lymphocyte activation antigen CD80, Activation B7-1 antigen, BB1, CTLA-4 counter-receptor B7.1, B7, CD antigen CD80, Activation B7-1 antigen, B lymphocyte activation antigen B7, B7,
Applications:	WB IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-2J8
Clonality:	Monoclonal Antibody
SwissProt ID:	P33681
Immunogen:	A synthesized peptide derived from human CD80
Gene ID:	941

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CD80

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 33 kDa, Observed MW: 60 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.