

Anti-TEF1/TEAD-1 [R02-2J3] Monoclonal Antibody

AGMB01626

Description

This Anti-TEF1/TEAD-1 [R02-2J3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01626
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	TEAD1, TCF13, TEF1, Transcriptional Enhancer Factor TEF 1, Transcriptional enhancer factor TEF-1, NTEF-1, Protein GT-IIC, TEA domain family member 1, TEAD-1, Transcription factor 13, TCF-13
Applications:	WB IHC-P IP
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R02-2J3
Clonality:	Monoclonal Antibody
SwissProt ID:	P28347
Immunogen:	A synthesized peptide derived from human TEAD1
Gene ID:	7003

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	TEAD1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 48 kDa, Observed MW: 52 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IP	1:50
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	