

Anti-CD105 [R05-3H7] Monoclonal Antibody

AGMB01704

Description

This Anti-CD105 [R05-3H7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01704
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ENG, END, ORW, HHT1, ORW1, CD105, FLJ41744, Endoglin, CD antigen CD105, CD105/Endoglin
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-3H7
Clonality:	Monoclonal Antibody
SwissProt ID:	P17813
Immunogen:	Recombinant protein of human CD105
Gene ID:	2022

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	ENG
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 71 kDa, Observed MW: 95 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
	ICC/IF	1:200-1:500
	FC	1:200-1:500
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.