

Anti-Spi1 [R02-1D7] Monoclonal Antibody

AGMB01718

Description

This Anti-Spi1 [R02-1D7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01718
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Dis1, PU.1, Dis-1, Sfpi1, Spi-1, Sfpi-1, Tcfpu1, Tfpu.1
Applications:	WB IHC-F IHC-P ICC/IF FC IP
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R02-1D7
Clonality:	Monoclonal Antibody
SwissProt ID:	P17433
Immunogen:	A synthetic peptide of mouse PU.1/Spi1
Gene ID:	20375
Gene Name:	Spi1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 31 kDa, Observed MW: 31 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:1000-1:5000
	ICC/IF	1:200-1:500
	FC	1:20-1:50
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.