

Anti-Interferon alpha/beta Receptor 1 [R09-7D6] Monoclonal Antibody

AGMB01723

Description

This Anti-Interferon alpha/beta Receptor 1 [R09-7D6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01723
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	IFNAR1, IFNAR, Interferon alpha/beta receptor 1, IFN-R-1, IFN-alpha/beta receptor 1, Cytokine receptor class-II member 1, Cytokine receptor family 2 member 1, CRF2-1, Type I interferon receptor 1
Applications:	WB FC IP
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R09-7D6
Clonality:	Monoclonal Antibody
SwissProt ID:	P17181
Immunogen:	A synthetic peptide of human Interferon alpha/beta receptor 1
Gene ID:	3454

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	IFNAR1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 64 kDa, Observed MW: 110-130 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IP	1:20

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.