

# Anti-PF4 [R09-7H9] Monoclonal Antibody

AGMB01787

## Description

---

This Anti-PF4 [R09-7H9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB01787
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	C-X-C motif chemokine 4, CXCL4, Iroplact, OncostatinA, PF4, Platelet factor 4, SCYB4, short form, Small inducible cytokine subfamily member 4, Chemokine (C X C motif) ligand 4, Chemokine (CXC motif) ligand 4,
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b> <b>IP</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human, Rat
<b>Clone:</b>	R09-7H9
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P02776
<b>Immunogen:</b>	A synthesized peptide derived from human PF4
<b>Gene ID:</b>	5196

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** PF4

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 11 kDa, Observed MW: 11 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IP	1:50

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.