

Anti-CD8 alpha [R07-6R-8] Monoclonal Antibody

AGMB01808

Description

This Anti-CD8 alpha [R07-6R-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|-----------------------|--|
| SKU: | AGMB01808 |
| Contents: | 50µl, 100µl Bradford Reagent: 1 vial (2ml) |
| Synonyms: | CD8A, MAL, T-cell surface glycoprotein CD8 alpha chain, T-lymphocyte differentiation antigen T8/Leu-2, CD antigen CD8a |
| Applications: | IHC-P FC |
| Research Area: | Immunology |
| Form: | Liquid |

Antibody Data

| | |
|----------------------|---------------------|
| Reactivity: | Human |
| Clone: | R07-6R-8 |
| Clonality: | Monoclonal Antibody |
| SwissProt ID: | P01732 |
| Immunogen: | Peptide |
| Gene ID: | 925 |

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CD8A
Host Species: Rabbit
Isotype: IgG
Purification: Protein A
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|---------------------------------|--------------------|--------------------------------|
| | IHC-P | 1:50-1:100 |
| | FC | 1:50 |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.