

Anti-TIMP1 [R03-9H5] Monoclonal Antibody

AGMB01828

Description

This Anti-TIMP1 [R03-9H5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01828
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Metalloproteinase inhibitor 1, Erythroid-potentiating activity (EPA), Fibroblast collagenase inhibitor (Collagenase inhibitor), Tissue inhibitor of metalloproteinases 1 (TIMP-1), TIMP1, CLGI, TIMP
Applications:	WB IHC-P ICC/IF
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-9H5
Clonality:	Monoclonal Antibody
SwissProt ID:	P01033
Immunogen:	A synthesized peptide derived from human TIMP1
Gene ID:	7076

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: TIMP1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 23 kDa, Observed MW: 26 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:100-1:200
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.