

Anti-Lactate Dehydrogenase A [R06-6J2] Monoclonal Antibody

AGMB01854

Description

This Anti-Lactate Dehydrogenase A [R06-6J2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01854
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	LDHA, LDH-M, LDH1, PIG19, GSD11, TRG5, Proliferation inducing gene 19
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-6J2
Clonality:	Monoclonal Antibody
SwissProt ID:	P00338
Immunogen:	A synthetic peptide of human Lactate Dehydrogenase
Gene ID:	3939
Gene Name:	LDHA

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 37 kDa, Observed MW: 37 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 1xPBS(pH 7.4), 150mM NaCl, 50% Glycerol, 0.02% Sodium azide and 0.05% BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:5000-1:150000
	IHC-P	1:200-1:2000
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	