

# Anti-NCR1 [R02-8W-5] Monoclonal Antibody

## AGMB01885

### Description

---

This Anti-NCR1 [R02-8W-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |                                                                                |
|-----------------------|--------------------------------------------------------------------------------|
| <b>SKU:</b>           | AGMB01885                                                                      |
| <b>Contents:</b>      | 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                  |
| <b>Synonyms:</b>      | natural cytotoxicity triggering receptor 1, LY94, CD335, NKP46, NK-p46, , NCR1 |
| <b>Applications:</b>  | <b>IHC-P</b>                                                                   |
| <b>Research Area:</b> | Immunology                                                                     |
| <b>Form:</b>          | Liquid                                                                         |

### Antibody Data

---

|                      |                                   |
|----------------------|-----------------------------------|
| <b>Reactivity:</b>   | Human                             |
| <b>Clone:</b>        | R02-8W-5                          |
| <b>Clonality:</b>    | Monoclonal Antibody               |
| <b>SwissProt ID:</b> | O76036                            |
| <b>Immunogen:</b>    | Recombinant protein of human NCR1 |
| <b>Gene ID:</b>      | 9437                              |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** NCR1

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | IHC-P              | 1:500-1:1000                   |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.