

Anti-KDM1/LSD1 [R08-3J5] Monoclonal Antibody

AGMB01905

Description

This Anti-KDM1/LSD1 [R08-3J5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB01905
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	KDM1A, AOF2, KDM1, KIAA0601, LSD1, Lysine-specific histone demethylase 1A, BRAF35-HDAC complex protein BHC110, Flavin-containing amine oxidase domain-containing protein 2
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R08-3J5
Clonality:	Monoclonal Antibody
SwissProt ID:	O60341
Immunogen:	A synthetic peptide of human KDM1/LSD1
Gene ID:	23028

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: KDM1A
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 93 kDa, Observed MW: 110 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC	1:50-1:200
	ICC/IF	1:100-1:200
	FC	1:20
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.