

Anti-CD79a [R02-7E3] Monoclonal Antibody

AGMB02047

Description

This Anti-CD79a [R02-7E3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02047
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD79A, IGA, MB1, B-cell antigen receptor complex-associated protein alpha chain, Ig-alpha, MB-1 membrane glycoprotein, Membrane-bound immunoglobulin-associated protein, Surface IgM-associated protein, CD antigen CD79a, B lymphocyte-specific MB1 protein, B-cell antigen receptor complex-associated protein alpha chain, CD 79a,
Applications:	WB IHC-F IHC-P ICC/IF
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R02-7E3
Clonality:	Monoclonal Antibody
SwissProt ID:	P11912
Immunogen:	A synthetic peptide of human CD79a

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 973

Gene Name: CD79A

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 25 kDa, Observed MW: 45-55 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:2000-1:5000
	IHC-P	1:100-1:200
	ICC/IF	1:200-1:500

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.