

Anti-CLASP1 [R01-2A2] Monoclonal Antibody

AGMB02064

Description

This Anti-CLASP1 [R01-2A2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU: AGMB02064

Contents: 20µl, 50µl, 100µl
Bradford Reagent: 1 vial (2ml)

Synonyms: CLASP1, KIAA0622, MAST1, CLIP-associating protein 1, Cytoplasmic linker-associated protein 1, Multiple asters homolog 1, Protein Orbit homolog 1, hOrbit1

Applications: **WB** **IHC-P** **ICC/IF** **FC**

Research Area: Cell Biology

Form: Liquid

Antibody Data

Reactivity: Human, Mouse, Rat

Clone: R01-2A2

Clonality: Monoclonal Antibody

SwissProt ID: Q7Z460

Immunogen: A synthetic peptide of human CLASP1

Gene ID: 23332

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CLASP1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 169 kDa, Observed MW: 160 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:200-1:500
	ICC/IF	1:50-1:100
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.