

Anti-Clusterin alpha Chain [R07-7A7] Monoclonal Antibody

AGMB02067

Description

This Anti-Clusterin alpha Chain [R07-7A7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02067
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Cli, ApoJ, Sgp2, SP-40, Sgp-2, Sugp-2, A1893575, D14Ucla3
Applications:	WB IHC-P IP
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Mouse, Rat
Clone:	R07-7A7
Clonality:	Monoclonal Antibody
SwissProt ID:	Q06890
Immunogen:	Recombinant protein of mouse Clusterin alpha chain
Gene ID:	12759
Gene Name:	Clu

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 52 kDa, Observed MW: 42 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IP	1:20
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.