

Anti-CRMP3 [R01-6G1] Monoclonal Antibody

AGMB02074

Description

This Anti-CRMP3 [R01-6G1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02074
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	DPYSL4, CRMP3, ULIP4, Dihydropyrimidinase-related protein 4, DRP-4, Collapsin response mediator protein 3, CRMP-3, UNC33-like phosphoprotein 4, ULIP-4
Applications:	WB
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R01-6G1
Clonality:	Monoclonal Antibody
SwissProt ID:	O14531
Immunogen:	A synthetic peptide of human CRMP3
Gene ID:	10570

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: DPYSL4

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 62 kDa, Observed MW: 62 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:2000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.