

Anti-Cytochrome C Oxidase subunit VIc/COX6C [R06-1D9] Monoclonal Antibody

AGMB02082

Description

This Anti-Cytochrome C Oxidase subunit VIc/COX6C [R06-1D9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02082
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	cytochrome c oxidase subunit 6C, COX6C, COX 6C, Cox6c, COX6C_HUMAN,
Applications:	WB IHC-P
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-1D9
Clonality:	Monoclonal Antibody
SwissProt ID:	P09669
Immunogen:	A synthetic peptide of human Cytochrome C Oxidase subunit Vic
Gene ID:	1345

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: COX6C

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 9 kDa, Observed MW: 12 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:200-1:500
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	