

Anti-Folate Binding Protein/FBP [R02-3G2] Monoclonal Antibody

AGMB02129

Description

This Anti-Folate Binding Protein/FBP [R02-3G2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02129
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FOLR1, FOLR, Folate receptor alpha, FR-alpha, Adult folate-binding protein, FBP, Folate receptor 1, Folate receptor, adult, KB cells FBP, Ovarian tumor-associated antigen MOv18
Applications:	WB IHC-P
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R02-3G2
Clonality:	Monoclonal Antibody
SwissProt ID:	P15328
Immunogen:	A synthetic peptide of human Folate Binding Protein
Gene ID:	2348

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: FOLR1

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 30 kDa, Observed MW: 40 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:200-1:1000
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	