

Anti-HLA A [R05-5H7] Monoclonal Antibody

AGMB02165

Description

This Anti-HLA A [R05-5H7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02165
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Aw-68, HLA class I histocompatibility antigen, A-28 alpha chain, MHC class I antigen A*68, HLA-A, MHC class I antigen HLA A heavy chain, ntigen presenting molecule, HLA class I histocompatibility antigen, A 1 alpha chain, HLAA,
Applications:	WB IHC-F IHC-P ICC/IF FC IP
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-5H7
Clonality:	Monoclonal Antibody
SwissProt ID:	P04439
Immunogen:	A synthetic peptide of human HLA A

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	3105
Gene Name:	HLA-A
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 41 kDa, Observed MW: 41 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:500-1:1000
IHC	1:50-1:100
ICC/IF	1:50-1:100
FC	1:10-1:100
IP	1:20

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.