

Anti-Ubiquitin K48 [R06-2H7] Monoclonal Antibody

AGMB02200

Description

This Anti-Ubiquitin K48 [R06-2H7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02200
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Ubiquitin K48, FLJ25987, MGC8385, ubiquitin B, Ubiquitin, UBCEP1, UBCEP2, RPS27A
Applications:	WB IHC-P ICC/IF FC
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-2H7
Clonality:	Monoclonal Antibody
SwissProt ID:	P0CG47
Immunogen:	A synthetic peptide of human K-48 Linkage Specific Ubiquitin
Gene ID:	7314

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: UBB

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:100-1:200
	ICC/IF	1:200-1:500
	FC	1:200-1:2000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.