

# Anti-LAMP2 [R03-2J8] Monoclonal Antibody

## AGMB02211

### Description

---

This Anti-LAMP2 [R03-2J8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB02211
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	LAMP2, Lysosome-associated membrane glycoprotein 2, LAMP-2, Lysosome-associated membrane protein 2, CD107 antigen-like family member B, CD107b, CD107 antigen-like family member B, CD107b, LAMP 2,
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b> <b>FC</b> <b>IP</b>
<b>Research Area:</b>	Tags & Cell Markers
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	R03-2J8
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P13473
<b>Immunogen:</b>	A synthetic peptide of human LAMP2A

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 3920  
**Gene Name:** LAMP2  
**Host Species:** Rabbit  
**Isotype:** IgG  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 45 kDa, Observed MW: 110 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

**Antibody Dilution Ratio:**

Application	Antibody Dilution Ratio
WB	1:500-1:1000
IHC	1:20-1:100
FC	1:100-1:200
IP	1:20

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.