

Anti-liver FABP [R08-7J1] Monoclonal Antibody

AGMB02215

Description

This Anti-liver FABP [R08-7J1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02215
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FABP 1, FABP1, FABP1, FABPL, FABPL_HUMAN, Fatty Acid Binding Protein 1, Fatty acid binding protein 1 liver, Fatty Acid Binding Protein, Fatty acid-binding protein 1, Fatty acid-binding protein, Fatty acid-binding protein liver, L FABP, L-FABP, liver, Liver-type fatty acid-binding protein.
Applications:	WB IHC-F IHC-P ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R08-7J1
Clonality:	Monoclonal Antibody
SwissProt ID:	P07148
Immunogen:	Recombinant protein of human liver FABP

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 2168
Gene Name: FABP1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 14 kDa, Observed MW: 14 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:500-1:1000
IHC-P	1:500-1:2000
ICC/IF	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.