

Anti-CSF1R/CD115 Receptor [R06-5H8] Monoclonal Antibody

AGMB02221

Description

This Anti-CSF1R/CD115 Receptor [R06-5H8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02221
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Macrophage colony-stimulating factor 1 receptor, CSF-1 receptor (EC:2.7.10.1), CSF-1-R, CSF-1R, M-CSF-R, Proto-oncogene c-Fms, CD115, Csf1r, Csfmr, Fms
Applications:	WB ICC/IF
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Mouse, Rat
Clone:	R06-5H8
Clonality:	Monoclonal Antibody
SwissProt ID:	P09581
Immunogen:	Recombinant protein of mouse MCSF Receptor
Gene ID:	12978

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: Csf1R

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 109 kDa, Observed MW: 109 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IF	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.