

Anti-NUDT19 [R05-9A5] Monoclonal Antibody

AGMB02263

Description

This Anti-NUDT19 [R05-9A5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02263
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Mitochondrial,Nudix motif 19,Nucleoside diphosphate-linked moiety X motif 19,Acyl-coenzyme A diphosphatase NUDT19,nudix (nucleoside diphosphate linked moiety X)-type motif 19,NUD19_HUMAN,RP2
Applications:	WB IHC-F IHC-P ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-9A5
Clonality:	Monoclonal Antibody
SwissProt ID:	A8MXV4
Immunogen:	Recombinant protein of human NUDT19
Gene ID:	390916

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NUDT19
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 42 kDa, Observed MW: 42 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:50-1:100
	ICC/IF	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.