

Anti-Nup107 [R06-9K7] Monoclonal Antibody

AGMB02265

Description

This Anti-Nup107 [R06-9K7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02265
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ODG6, NUP84, NPHS11, ODG6, GAMOS7, 107 kDa nucleoporin, NU107_HUMAN, Nuclear pore complex protein Nup107,
Applications:	
Research Area:	Tags & Cell Markers
Form:	Liquid

Antibody Data

Reactivity:	Human, Rat
Clone:	R06-9K7
Clonality:	Monoclonal Antibody
SwissProt ID:	P57740
Immunogen:	A synthetic peptide of human Nup107
Gene ID:	57122

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NUP107
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 106 kDa, Observed MW: 106 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:2000-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.