

**Anti-HAUS7 [R05-9C8] Monoclonal Antibody**  
**AGMB02398**

**Description**

---

This Anti-HAUS7 [R05-9C8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

**Product Information**

---

|                |                                                                                     |
|----------------|-------------------------------------------------------------------------------------|
| SKU:           | AGMB02398                                                                           |
| Contents:      | 20µl, 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml)                                 |
| Synonyms:      | UIP1, UCHL5IP                                                                       |
| Applications:  |  |
| Research Area: | Cell Biology                                                                        |
| Form:          | Liquid                                                                              |

**Antibody Data**

---

|               |                                      |
|---------------|--------------------------------------|
| Reactivity:   | Human                                |
| Clone:        | R05-9C8                              |
| Clonality:    | Monoclonal Antibody                  |
| SwissProt ID: | Q99871                               |
| Immunogen:    | Recombinant protein of human UCHL5IP |
| Gene ID:      | 55559                                |
| Gene Name:    | HAUS7                                |

**Manufacturers Statement**

This final kit system is assembled and quality-released by Assay Genie Limited.

|                          |                                            |
|--------------------------|--------------------------------------------|
| <b>Host Species:</b>     | Rabbit                                     |
| <b>Isotype:</b>          | IgG                                        |
| <b>Purification:</b>     | Affinity Purified                          |
| <b>Conjugated:</b>       | Unconjugated                               |
| <b>Modification:</b>     | Unmodified                                 |
| <b>Molecular Weight:</b> | Calculated MW: 40 kDa, Observed MW: 40 kDa |

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

| <b>Antibody Dilution Ratio:</b> | <b>Application</b> | <b>Antibody Dilution Ratio</b> |
|---------------------------------|--------------------|--------------------------------|
|                                 | WB                 | 1:2000-1:5000                  |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.