

Anti-CLPX [R05-2A6] Monoclonal Antibody

AGMB02461

Description

This Anti-CLPX [R05-2A6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02461
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	ATP-dependent Clp protease ATP-binding subunit clpX-like, mitochondrial, CLPX, ATP dependent Clp protease ATP binding subunit clpX like, mitochondrial, ATP-dependent Clp protease ATP-binding subunit clpX-like, mitochondrial, AU014732,
Applications:	WB ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R05-2A6
Clonality:	Monoclonal Antibody
SwissProt ID:	O76031
Immunogen:	A synthetic peptide of human CLPX

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 10845
Gene Name: CLPX
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 69 kDa, Observed MW: 69 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:1000-1:5000
ICC/IF	1:100-1:200

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.