

# Anti-Transcriptional Enhancer Factor TEF [R05-5P3] Monoclonal Antibody

## AGMB02648

### Description

---

This Anti-Transcriptional Enhancer Factor TEF [R05-5P3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB02648
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	DTEF 1, NTEF 1, TCF13, TEAD1, TEAD2, TEAD3, TEAD4, TEAD5, TEF1, TEF4, TEF5
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Epigenetics
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R05-5P3
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P28347/Q15561/Q15562/Q99594
<b>Immunogen:</b>	A synthesized peptide derived from human TEAD1 + TEAD2 + TEAD3 + TEAD4

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 7003/7004/8463/7005

**Gene Name:** TEAD1/2/3/4

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 48 kDa, Observed MW: 54 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.