

Anti-Guanylyl Cyclase beta 1 [R05-3C1] Monoclonal Antibody

AGMB02652

Description

This Anti-Guanylyl Cyclase beta 1 [R05-3C1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB02652
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	GCbeta1, GCS beta 1, GCS beta 3, GUC1B3, GUCSB3, GUCY1B1, Gucy1b1, GUCY1B3, SGC, GC-S-beta-1, GC-SB3, GCbeta1,
Applications:	WB IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R05-3C1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q02153
Immunogen:	A synthesized peptide derived from human Guanylyl Cyclase beta 1
Gene ID:	2983

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	GUCY1B1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 71 kDa, Observed MW: 71 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	